

# Han Li Lab PyMol Guide

University of California, Irvine

March 31, 2020

# Table of Contents

- ① Object Selection
- ② Structure Manipulation
- ③ Publication Quality Images
- ④ Rotation Movies

# Selection Shortcuts

<code>sel &lt;variable_name&gt;, &lt;object&gt;</code>	<code># selection template</code>
<code>remove solv</code>	<code># remove solvent</code>
<code>remove inorg</code>	<code># remove ions</code>
<code>remove not chain &lt;chain&gt;</code>	<code># remove other chains</code>
<code>sel &lt;ligand&gt;, org</code>	<code># select all ligands</code>
<code>sel &lt;protein&gt;, poly</code>	<code># select all protein atoms</code>
<code>sel &lt;target&gt;, rep sticks</code>	<code># select by representation</code>
<code>sel &lt;residue&gt;, resi &lt;residue number&gt;</code>	<code># select residue by number</code>
<code>sel &lt;residue&gt;, resn &lt;residue name&gt;</code>	<code># select object by residue name</code>
<code>sel &lt;atom&gt;, name &lt;atom name&gt;</code>	<code># select atom by atom name</code>
<code>reset</code>	<code># reset camera origin</code>

## `# Examples`


`# selects residues 23, 54, 55, 56, 57, 99`

`sel target_residues, resi 23+54-57+99`

`# selects ligand with 3-letter code NAD from chain A`

`sel target_ligand, resn NAD and chain A`

# Displaying Binding Pocket Interactions



```
# all residues with any atom  
    around 3.5 angstrom  
sel pocket, br. org around 3.5  
show sticks, pocket  
color cyan, pocket  
util.cnc  
dist hbonds, org, poly, mode=2  
hide labels
```

# Alignment

```
# 1 to 1 alignment
```

```
fetch 1j49; fetch 4e5n
```

```
remove solv; remove not chain A; remove inorg
```

```
align 1j49, 4e5n # fits 1j49 onto 4e5n
```

```
reset
```

```
# multiple to 1 alignment
```

```
fetch 1j49; fetch 4e5n; fetch 6ih4
```

```
remove solv; remove not chain A; remove inorg
```

```
align to 4e5n # align all to 4e5n
```

```
reset
```

# Ligand Transfer

```
fetch 1qi1; fetch 1qi6           # 1qi1 gapn with nadp, 1qi6 apo gapn
remove not chain A
remove solv
remove inorg
align 1qi1, 1qi6                 # align binding pockets
extract nadp, resn NAP
sel 1qi6_nadp, 1qi6 + nadp
save 1qi6_nadp.pdb, 1qi6_nadp    # or save with File -> Export Molecule
# type pwd to see the default directory the file is saved to
```

# Mutagenesis

Follow instructions from the PyMol website

<https://pymolwiki.org/index.php/Mutagenesis>

- ① Wizard → Mutagenesis → Protein
- ② Select target residue by clicking
- ③ Select mutation residue in Mutagenesis box
- ④ Examine rotamers with arrows on bottom right
- ⑤ Apply to save



# General Settings

```
bg white                                # change background color
set ray_opaque_background, 1           # make background solid
set ray_shadow, 0                      # turn off shadows
set antialias, 2                      # smooth edges
set cartoon_fancy_helices, 1          # adds ridges to helix edges
util.cnc                              # colors oxygen and nitrogen atoms
ray <width>                            # ray and scale resolution

# Recommended settings
set ray_trace_mode, 3                  # optional neon colors and outline
color white, poly
color green, org
color cyan, poly and rep sticks
```




## Showing Backbone Atoms




```
fetch 1j49  
[...] # clean up  
dist hbonds, org, resi 234, mode=2  
hide labels  
set cartoon_side_chain_helper, 0, resi 234
```

# Glucose Dehydrogenase




```
load gludh_80_nmn.pdb
bg white
set ray_opaque_background, 1
set cartoon_fancy_helices, 1
set ray_shadow, 0
color white, poly
color orange, org
show cart
show surf
set transparency, 0.7
sel mutations, i. 93+39+195+17 and
    name ca
show spheres, mutations
color deeppurple, rep spheres
set sphere_scale, 0.5
util.cnc org
```

# Lactate Dehydrogenase



```
fetch 1j49
bg white; color white, poly
remove chain A; remove inorg;
    remove solv
set ray_opaque_background, 1
set ray_trace_mode, 3
set cartoon_transparency, 0.3
color green, org
show lines, br. org around 3.5
color silver, rep lines
sel spec_loop, resi 176-178
show sticks, spec_loop
color cyan, resi 176
color violet, resi 177
color lightorange, resi 178
dist hbonds, resi 176, all, mode=2
hide labels; show lines, resi 154
set cartoon_side_chain_helper, 0,
    resi 154+177+178
util.cnc
```

# Nitrogenase Iron Protein



```
fetch 6nzj
bg white
color lightblue, chain A
color palegreen, chain B
util.cnc
remove resn S04; remove solv
set ray_opaque_background, 1
# outline with regular colors
set ray_trace_mode, 1
set cartoon_transparency, 0.7
show sticks, inorg
set sphere_scale, 0.4
set cartoon_transparency, 0.5
show sticks, br. inorg around 3.5
hide sticks, resi 96
```

# Rotation Movies

Common options under Movie → Program

GluDH

```
movie.add_roll(16.0, axis='y', start=1)
```

File -> Export Movie

<https://youtu.be/eaRJjZ0OLRs>